

A dark, close-up photograph of a dental procedure. A handpiece is visible on the left, and a hand is holding a dental instrument on the right, positioned near a patient's teeth. The background is dark and out of focus.

DYNAMIC
DENTAL EQUIPMENT

DIVERSE BURR REPLACEMENT

A close-up, dimly lit photograph of a dental procedure. A dental professional's hand is visible, using a specialized tool to adjust a component on a patient's head. The patient's head is tilted back, and the dental equipment is metallic and complex. The overall tone is professional and clinical.

YOU WILL NEED:

AVAILABLE FROM YOUR DYNAMIC DISTRIBUTOR:

- DYNAMIC HEAD LOCKING TOOL
- 21MM SPANNER
- DYNAMIC NEW REPLACEMENT BURR

- UNSCREW THE NUT BEHIND THE HEAD IN AN ANTI-CLOCKWISE DIRECTION. IF STIFF USE 21MM SPANNER TO LOOSEN.



▪ SCREW THE HEAD LOCKING
TOOL ON TO THE BACK OF THE
HEAD CASING



■ UNSCREW BURR ANTI-
CLOCKWISE, AND REMOVE






- SCREW IN NEW BURR
CLOCKWISE UNTIL TIGHT,
BUT DO NOT OVER TIGHTEN


- UNSCREW THE HEAD LOCKING TOOL ON THE BACK OF THE HEAD CASING





- APPLY A LITTLE FOOD LUBE GREASE AROUND THE MALE X-DRIVE AND THREAD AS SHOWN

- ATTACH HEAD BACK ONTO THE SHAFT BY TIGHTENING THE NUT ON THE TOP OF THE SHAFT IN A CLOCKWISE DIRECTION. TAKING CARE NOT TO CROSS THREAD.
- **IMPORTANT ONLY HAND TIGHTEN**



NOTE: SEAL PRESENT